

Project Details

|           |  |
|-----------|--|
| Name:     | taenia_solium_prjna170813_wbps19_cds_transcripts |
| Creation: | Wed Sep 11 2024 11:13:34                         |
| Version:  | OmicsBox 3.3.2                                   |

Analysis Steps

|                               |                          |
|-------------------------------|--------------------------|
| Blast Description Annotator   |                          |
| Modification:                 | Fri Sep 13 2024 03:12:56 |
| Version:                      | OmicsBox 3.3.2           |
| Parameters:                   |                          |
| Recover original description: | false                    |
| Comment:                      |                          |

|                               |                          |
|-------------------------------|--------------------------|
| Blast Description Annotator   |                          |
| Modification:                 | Fri Sep 13 2024 02:52:57 |
| Version:                      | OmicsBox 3.3.2           |
| Parameters:                   |                          |
| Recover original description: | true                     |
| Comment:                      |                          |

|                                      |                          |
|--------------------------------------|--------------------------|
| Merge InterProScan GOs to Annotation |                          |
| Modification:                        | Fri Sep 13 2024 02:27:14 |
| Version:                             | OmicsBox 3.3.2           |
| Comment:                             |                          |

|                                      |                          |
|--------------------------------------|--------------------------|
| Annotate GOs from Blast Descriptions |                          |
| Modification:                        | Fri Sep 13 2024 02:26:20 |
| Version:                             | OmicsBox 3.3.2           |
| Parameters:                          |                          |
| Validate Annotations:                | true                     |
| Minimum Similarity:                  | 85.0                     |
| Comment:                             |                          |

|                     |                          |
|---------------------|--------------------------|
| Validate Annotation |                          |
| Modification:       | Fri Sep 13 2024 02:25:38 |
| Elapsed Time:       | 10s                      |
| Version:            | OmicsBox 3.3.2           |
| Comment:            |                          |

|                       |                          |
|-----------------------|--------------------------|
| GO-EnzymeCode Mapping |                          |
| Modification:         | Fri Sep 13 2024 02:24:09 |
| Version:              | OmicsBox 3.3.2           |
| Comment:              |                          |

**Merge EggNOG 5 GO Annotations**

|                                 |   |
|---------------------------------|---|
| Modification:                   | Fri Sep 13 2024 02:20:11  |
| Elapsed Time:                   | 12s   |
| Version:                        | OmicsBox 3.3.2  |
| Parameters:                     |   |
| Seed Ortholog E-Value Filter:   | 1E-5  |
| Seed Ortholog Bit-Score Filter: | 60.0  |
| EggNOG Annotations:             | C:<br>\Users\dcc_9\OneDrive\Escritorio\Anno\EggNOG\taenia_solium_prjna170813_wbps19_cds_transcripts_fa_eggbox |
| Comment:                        |   |

**Annotation**

|                          |                          |
|--------------------------|--------------------------|
| Modification:            | Fri Sep 13 2024 02:11:24 |
| Version:                 | OmicsBox 3.3.2           |
| Parameters:              |                          |
| E-Value-Hit-Filter:      | 1.0E-6                   |
| Annotation CutOff:       | 55                       |
| NAS:                     | 0.8                      |
| IPI:                     | 1.0                      |
| NR:                      | 0.0                      |
| IDA:                     | 1.0                      |
| TAS:                     | 0.9                      |
| IBA:                     | 0.8                      |
| IBD:                     | 0.8                      |
| Hit Filter:              | 500                      |
| HSP-Hit Coverage CutOff: | 0                        |
| ISA:                     | 0.8                      |
| IC:                      | 0.9                      |
| ISM:                     | 0.8                      |
| ISO:                     | 0.8                      |
| IGC:                     | 0.7                      |
| IEA:                     | 0.7                      |
| GO Weight:               | 5                        |
| Taxon:                   |                          |
| ISS:                     | 0.8                      |
| IGI:                     | 1.0                      |
| IMP:                     | 1.0                      |
| Only hits with GOs:      | false                    |
| IKR:                     | 0.8                      |
| RCA:                     | 0.8                      |
| IEP:                     | 1.0                      |
| ND:                      | 0.5                      |
| IRD:                     | 0.7                      |
| EXP:                     | 1.0                      |
| Comment:                 |                          |

| GO-EnzymeCode Mapping |                          |
|-----------------------|--------------------------|
| Modification:         | Fri Sep 13 2024 02:11:02 |
| Version:              | OmicsBox 3.3.2           |
| Comment:              |                          |

| Merge EggNOG 5 GO Annotations   |   |
|---------------------------------|---|
| Modification:                   | Fri Sep 13 2024 02:01:05  |
| Elapsed Time:                   | 13s   |
| Version:                        | OmicsBox 3.3.2  |
| Parameters:                     |   |
| Seed Ortholog E-Value Filter:   | 1E-5  |
| Seed Ortholog Bit-Score Filter: | 60.0  |
| EggNOG Annotations:             | C:\Users\dcc_9\OneDrive\Escritorio\Anno\EggNOG\taenia_solium_prjna170813_wbps19_cds_transcripts_fa_eggbox |
| Comment:                        |   |

| Annotation               |                          |
|--------------------------|--------------------------|
| Modification:            | Fri Sep 13 2024 02:00:08 |
| Elapsed Time:            | 10s                      |
| Version:                 | OmicsBox 3.3.2           |
| Parameters:              |                          |
| E-Value-Hit-Filter:      | 1.0E-6                   |
| Annotation CutOff:       | 55                       |
| NAS:                     | 0.8                      |
| IPI:                     | 1.0                      |
| NR:                      | 0.0                      |
| IDA:                     | 1.0                      |
| TAS:                     | 0.9                      |
| IBA:                     | 0.8                      |
| IBD:                     | 0.8                      |
| Hit Filter:              | 500                      |
| HSP-Hit Coverage CutOff: | 0                        |
| ISA:                     | 0.8                      |
| IC:                      | 0.9                      |
| ISM:                     | 0.8                      |
| ISO:                     | 0.8                      |
| IGC:                     | 0.7                      |
| IEA:                     | 0.7                      |
| GO Weight:               | 5                        |
| Taxon:                   | 6157 Platyhelminthes     |
| ISS:                     | 0.8                      |
| IGI:                     | 1.0                      |
| IMP:                     | 1.0                      |
| Only hits with GOs:      | false                    |
| IKR:                     | 0.8                      |
| RCA:                     | 0.8                      |
| IEP:                     | 1.0                      |
| ND:                      | 0.5                      |
| IRD:                     | 0.7                      |
| EXP:                     | 1.0                      |
| <a href="#">Comment:</a> |                          |

| Mapping                      |                          |
|------------------------------|--------------------------|
| Modification:                | Fri Sep 13 2024 01:54:38 |
| Elapsed Time:                | 3m17s                    |
| Version:                     | OmicsBox 3.3.2           |
| Parameters:                  |                          |
| Use latest database version: | true                     |
| Goa Version:                 | 2024_05                  |
| <a href="#">Comment:</a>     |                          |

| InterProScan             |   |
|--------------------------|---|
| Modification:            | Thu Sep 12 2024 11:55:07                            |
| Elapsed Time:            | 36h28m10s   |
| Version:                 | OmicsBox 3.3.2                                      |
| Parameters:              |   |
| Coils:                   | true  |
| Choose Folder:           | C:\Users\dcc_9\OneDrive\Escritorio\Anno\IPS-Results |
| CDD:                     | true  |
| Email:                   | daccar98@gmail.com                                  |
| Superfamily:             | true  |
| Panther:                 | true  |
| FunFam:                  | true  |
| PIRSF:                   | true  |
| ProSitePatterns:         | true  |
| Phobius:                 | true  |
| PfamA:                   | true  |
| NCBIfam:                 | true  |
| Gene3D:                  | true  |
| TMHMM:                   | true  |
| SFLD:                    | true  |
| AntiFam:                 | true  |
| PRINTS:                  | true  |
| ProSiteProfiles:         | true  |
| PIRSR:                   | true  |
| HAMAP:                   | true  |
| Export Formats:          | tsv,gff,xml   |
| SignalP:                 | true  |
| SMART:                   | true  |
| <a href="#">Comment:</a> |   |

| Blast Description Annotator   |                          |
|-------------------------------|--------------------------|
| Modification:                 | Wed Sep 11 2024 12:24:38 |
| Version:                      | OmicsBox 3.3.2           |
| Parameters:                   |                          |
| Recover original description: | true                     |
| <a href="#">Comment:</a>      |                          |

| Blast Description Annotator   |                          |
|-------------------------------|--------------------------|
| Modification:                 | Wed Sep 11 2024 12:24:21 |
| Version:                      | OmicsBox 3.3.2           |
| Parameters:                   |                          |
| Recover original description: | false                    |
| <a href="#">Comment:</a>      |                          |

| Diamond Blast                      |  |
|------------------------------------|--|
| Modification:                      | Wed Sep 11 2024 12:21:01   |
| Elapsed Time:                      | 58m46s   |
| Version:                           | OmicsBox 3.3.2   |
| Parameters:                        |  |
| Optimize for Long Query Seqs:      | true   |
| HSP-Hit Coverage:                  | 0  |
| HSP Length Cutoff:                 | 33   |
| Filter option:                     | Include only the selected taxonomies in the search   |
| Database:                          | NR (2024-07-11)  |
| Taxonomy Filter:                   | 6157 Platyhelminthes, 6178 Trematoda, 6199 Cestoda, 6231 Nematoda, 6239 Caenorhabditis elegans |
| Sensitivity Mode:                  | More Sensitive   |
| Blast Expectation Value (e-Value): | 1.0E-5   |
| Blast Mode:                        | blastx   |
| Comment:                           |  |

| Diamond Blast                      |  |
|------------------------------------|--|
| Modification:                      | Wed Sep 11 2024 12:21:01   |
| Elapsed Time:                      | 58m46s   |
| Version:                           | OmicsBox 3.3.2   |
| Parameters:                        |  |
| Optimize for Long Query Seqs:      | true   |
| HSP-Hit Coverage:                  | 0  |
| HSP Length Cutoff:                 | 33   |
| Filter option:                     | Include only the selected taxonomies in the search   |
| Database:                          | NR (2024-07-11)  |
| Taxonomy Filter:                   | 6157 Platyhelminthes, 6178 Trematoda, 6199 Cestoda, 6231 Nematoda, 6239 Caenorhabditis elegans |
| Sensitivity Mode:                  | More Sensitive   |
| Blast Expectation Value (e-Value): | 1.0E-5   |
| Blast Mode:                        | blastx   |
| Comment:                           |  |

| InterProScan             |   |
|--------------------------|---|
| Modification:            | Wed Sep 11 2024 11:26:56                            |
| Version:                 | OmicsBox 3.3.2                                      |
| Parameters:              |   |
| Coils:                   | true  |
| Choose Folder:           | C:\Users\dcc_9\OneDrive\Escritorio\Anno\IPS-Results |
| CDD:                     | true  |
| Email:                   | daccar98@gmail.com                                  |
| Superfamily:             | true  |
| Panther:                 | true  |
| FunFam:                  | true  |
| PIRSF:                   | true  |
| ProSitePatterns:         | true  |
| Phobius:                 | true  |
| PfamA:                   | true  |
| NCBIfam:                 | true  |
| Gene3D:                  | true  |
| TMHMM:                   | true  |
| SFLD:                    | true  |
| AntiFam:                 | true  |
| PRINTS:                  | true  |
| ProSiteProfiles:         | true  |
| PIRSR:                   | true  |
| HAMAP:                   | true  |
| Export Formats:          | tsv,gff,xml   |
| SignalP:                 | true  |
| SMART:                   | true  |
| <a href="#">Comment:</a> |   |

| Load Sequences                |   |
|-------------------------------|---|
| Modification:                 | Wed Sep 11 2024 11:13:38  |
| Version:                      | OmicsBox 3.3.2  |
| Parameters:                   |   |
| Fasta File:                   | C:\Users\dcc_9\OneDrive\Escritorio\Anno\Transcripts\taenia_solium.PRJNA170813.WBPS19.CDS_transcripts.fa |
| How to Handle Duplicated IDs: | Replace   |
| Add the Description:          | true  |
| <a href="#">Comment:</a>      |   |

| Load Sequences                |   |
|-------------------------------|---|
| Modification:                 | Wed Sep 11 2024 11:13:38  |
| Version:                      | OmicsBox 3.3.2  |
| Parameters:                   |   |
| Fasta File:                   | C:\Users\dcc_9\OneDrive\Escritorio\Anno\Transcripts\taenia_solium.PRJNA170813.WBPS19.CDS_transcripts.fa |
| How to Handle Duplicated IDs: | Replace   |
| Add the Description:          | true  |
| <a href="#">Comment:</a>      |   |